



SAMPLE
TRANSPORT
PACKAGING









The PathoPak range provides strong, durable and reusable, rigid secondary containers, compliant with transport regulations for UN Class 6.2 Infectious Substances Category A and B. Available in three sizes 1, 2 and 3 Litres.

THIS PRODUCT CAN BE USED FOR

-  Infectious samples
-  Diagnostic samples
-  Clinical trials
-  Screening programmes
-  Eradication schemes

FEATURES:

-  Complete Category A and Category B solutions include 95kPa tested secondary container, absorbent, bubble wrap, labels, and UN markings
-  Reusable
-  Ergonomically engineered bottle design for safe manual handling including 1.5 turn lid and non-tapered aperture for easy extraction of samples
-  4GU certified for the use of any primary
-  Suitable for cleaning by autoclave
-  Full refurbishment services available





PATHOPAK RANGE

Product Code	Product Name	External Dimensions	Payload Dimensions (Bottle opening diameter x internal height)	Payload Volume	Mass	Volumetric Weight L/6
PP011	DGP PathoPak 1L	112 x 110 x 190 mm 4.4 x 4.3 x 7.5 in	90 x 170 mm 3.5 x 6.7 in	1 L	0.24 kg 0.53 lbs	0.39 kgvol ^c 0.86 lbsvol ^c
PP012	DGP PathoPak 2L	153 x 156 x 193 mm 6.0 x 6.1 x 7.6 in	125 x 170 mm 4.9 x 6.7 in	2 L	0.45 kg 0.99 lbs	0.77 kgvol ^c 1.69 lbsvol ^c
PP013	DGP PathoPak 3L	160 x 160 x 280 mm 6.3 x 6.3 x 11.0 in	125 x 255 mm 4.9 x 10.0 in	3 L	0.6 kg 1.32 lbs	1.19 kgvol ^c 2.63 lbsvol ^c